
ENDOSKOP 401-B

User Manual

Čeština

Română

Polski

Slovenčina

Deutsch

Slovenščina

Magyar

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
English

Sunnysoft s.r.o.

401-B

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Obsah

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
User guide

Contents

1	Product description	3
2	Safety instructions	3
3	Hygiene and cleaning	4
4	Disposal	4
5	Controls and parts	4
6	Charging	5
7	Switching on and off	5
8	Connecting to a mobile phone	5
9	What the app can do	6
10	Connecting to a computer	6
11	Troubleshooting	7
12	Technical data	7
13	Conformity	8

1. Product description

Thank you for choosing the 401-B oral endoscope. It is a small inspection camera with its own lighting, intended for viewing places you cannot see into, chiefly the oral cavity.

The camera connects wirelessly to a mobile phone or by cable to a computer, and you watch the picture on their screen. The device itself has no display.

The device is for viewing and for taking photographs and video. It does not replace an examination by a dentist and is not intended for making a diagnosis.

Read this guide before you use the device for the first time, and keep it.

2. Safety instructions

Warning

Insert the probe into the mouth carefully and without force. Do not press it against teeth, gums or the mucous membrane.

- Do not use the device on injured or inflamed mucous membrane, or at the site of a recent procedure.
- If you see something in your mouth that worries you, consult a dentist. Do not treat yourself on the basis of the camera picture.
- Do not look directly into the lit LEDs from close up, and do not shine them into anyone's eyes.
- Only the probe head with the lens is waterproof. Do not immerse the body of the device, the buttons or the stand in water, and do not rinse them under running water.
- The device is not a toy. Give it to children only under adult supervision.
- Do not expose the device to temperatures outside the range given in the technical data, to direct sunlight or to lasting damp.
- Do not modify or repair the device yourself.

Warning

The battery is built in and cannot be removed. Do not dismantle the device, puncture it, throw it into fire or expose it to heat. Charge it only from an adapter of matching rating.

3. Hygiene and cleaning

- Clean the probe head before and after every use.
- Use a cotton swab moistened with water or alcohol and wipe the lens gently. Let it dry before you put the device away.
- Wipe the body of the device and the stand with a damp cloth only. Do not use solvents or abrasives: a scratched lens will no longer give a sharp picture.
- If several people use the device, clean the probe head between users.
- Do not boil the probe, do not sterilise it in an autoclave and do not wash it in a dishwasher.

4. Disposal


The device contains a built-in lithium battery that cannot be removed. It does not belong in household waste or in fire; hand it in whole at a collection point.

This symbol on the product or its packaging means that the product does not belong in household waste. Hand it in at a collection point for electrical and electronic equipment. Correct disposal helps prevent harm to health and to the environment.

5. Controls and parts

The numbering matches the illustration.

1. **Lens** with eight LEDs and a high-resolution camera.
2. **Photo and video:** a short press takes a photograph, a press of 3 seconds starts video recording.
3. **Zoom:** a short press cycles the magnification $1.0\times \rightarrow 1.5\times \rightarrow 2.0\times \rightarrow 1.0\times$.
4. **Lighting:** a short press switches the brightness between 100 %, 80 %, 50 % and off. A press of 3 seconds with the device switched on changes to computer-connection mode.
5. **WiFi indicator:** flashes until the connection is established; once connected it stays lit.
6. **Battery indicator:** flashes red when the battery is flat, lights red while charging and green when fully charged.
7. **Power button:** a press of 3 seconds switches the device on and off.

6. Charging

When the battery is running low, the battery indicator flashes red. Charge the device from an ordinary USB mains adapter of the rating given in the technical data.

While charging the indicator lights red, and green once fully charged. A full charge takes about 2 hours. Charge the device fully before first use.

7. Switching on and off

Press the power button and hold it for 3 seconds. The device switches on and starts broadcasting its own WiFi network.

The same press of 3 seconds switches the device off.

8. Connecting to a mobile phone

1. Install the ANESOK app. On an Android phone you will find it in the Google Play store, on an Apple phone in the App Store.
2. Switch the device on by pressing the power button for 3 seconds.

3. In the phone settings, open the list of WiFi networks.
4. Connect to the network whose name begins with ANESOK-401. The rest of the name differs from one unit to another.
5. Start the ANESOK app. The camera picture appears on the phone screen.

Note

While the phone is connected to the device's network it has no internet access. That is normal: the device creates its own network but does not connect it to anything else.

9. What the app can do

The app gives you the following functions:

- rotating the picture
- setting the resolution
- saved photographs and videos
- switching between colour and black-and-white
- the device's battery level
- comparing two shots side by side
- taking a photograph and recording video
- zooming in and out

10. Connecting to a computer

The device connects to a computer with a USB-C cable, not over WiFi. The procedure is the same in both systems apart from the program used to view the picture.

1. Connect the device to the computer with a USB-C cable and switch it on.
 2. Press the lighting button and hold it for 3 seconds. The WiFi indicator goes out, which means the device has changed to computer camera mode.
- **Windows:** start a program that displays the picture from a USB camera.
 - **macOS:** start Photo Booth, choose the USB camera from the Camera menu, and the picture appears.

Note

For Windows the manufacturer points to a program called Smart Camera. We do not repeat the download link given on the leaflet: it leads to an executable file hosted by a third party under a different brand name. Any program that supports a USB camera will display the picture.

11. Troubleshooting

Problem	What to do
The picture is blurred or unclear	Wipe the lens with a cotton swab moistened with water or alcohol.
The device will not switch on	Charge the battery.
The app has stopped responding	Close the apps running on the phone and switch the device on again.
The lighting does not come on	Restart the device.
The indicator is lit red	The battery is flat. Charge the device before use.
Everything is connected and the picture is running, but it does not match where the camera is pointing	Disconnect the phone from the device's network, restart the device, open the app again and establish the connection afresh.
The WiFi connection is unstable and drops after a few seconds	Restart the device and connect to its network again.

If the problems persist, contact your retailer.

12. Technical data

Model	401-B
Type	oral inspection camera with WiFi
Network	IEEE 802.11 b/g/n
Operating band	2.4 GHz
Antenna	built-in
Image resolution	1920 × 1080
Frame rate	25 fps
Image sensor	CMOS

Lighting	8 adjustable LEDs around the lens
Field of view	70°
Operating temperature	0 – 45 °C
Battery	lithium-ion, built-in, 3.7 V, 850 mAh
Runtime per charge	3 h
Power supply	5 V DC, 500 mA / 1 A
Ingress protection	IP67 (probe head with lens only)

Note

According to the manufacturer, the IP67 rating applies to the probe head with the lens only — the device as a whole is not waterproof.

13. Conformity

CE With the CE marking the manufacturer declares that the product meets the requirements of the European Union rules that apply to it.

The source names neither the manufacturer nor its address. The distributor's details are in the footer.

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